

Computer Requirements

MedRx **Minimum** Computer Specs:

- Windows® - PC compatible computer
- Intel™ i5, 2.0 GHz or better
- 4 GB RAM
- 20 GB free hard drive space
- Available 2.0 USB port
- Windows 10 or 11 Professional
- Compatible with 3.0 USB

MedRx **Recommended** Computer Specs:

- Windows® - PC Computer
- Intel™ i5 Quad Core or better
- 8 GB RAM or more
- Available USB 3.0
- Graphics Adapter with 2GB Dedicated Video Memory
- 50 GB or more free hard drive space
- High Speed Internet Connection
- Windows 10 or 11 Professional
- Compatible with 3.0 USB